



## *Certificate of Analysis*

*Catalog Number:* 1027  
*Lot Number:* 389AP10  
*Cell count:* 10 x 10<sup>6</sup> per vial  
*Viability:* 92%  
*CD3 positive:* 10.87  
*CD19 positive:* 5.9  
*CD14 positive:* 9.8  
*CD11c positive:* 65.0  
*CD16 positive:* 53.7  
*CD56 positive:* 81.3


### *Donor Information*

*Donor ID:* 091  
*ABO Type:* AB positive  
*Gender:* Male  
*Race:* African-American  
*Age:* 40

### *Notes:*

Donors are tested for the blood borne pathogens HIV-1 and 2, Hepatitis B, Hepatitis C and HTLV-1 and are negative. Cells should still be handled as if potentially infectious following biosafety level 2 procedures.

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE



**Astarte**  
BIOLOGICS

NK cells Positively Selected  
Catalog # 1027  
Lot # 389AP10

