



Certificate of Analysis

Catalog Number:	1021
Description:	Chondrocytes (Bovine) Primary Cells from 1 month old
Lot Number:	2700C09
Cell Count:	approximately 10×10^6 /vial
Viability:	80.0%
<u>Performance Data</u>	



Chondrocytes were thawed and cultured in 96 well plates at 10^6 cells/ml for 2-4 days. Culture supernatants were tested for proteoglycan using the Rheumera[®] GAG Assay Kit. See www.astartebio.com for more details.



Chondrocytes were thawed and cultured in 96 well plates at 10^6 cells/ml for 4 days. Type II collagen was measured using the Rheumera[®] Type II Collagen Capture ELISA Kit. For more experimental details see www.astartebio.com.