



Certificate of Analysis

Catalog Number: 9000
Description: Rat type II collagen immunized spleenocytes
Lot Number: 83MY07
Cell Count: 70 Million cells per vial
Viability: 80%

Donor Information

Species: Rat
Strain: LEWIS
Gender: Male
Notes: Rats immunized and boosted with native chick Type II collagen

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE

Astarte Biologics, LLC 866.513.1894 www.astartebio.com